

# MODELSCART

## CPR and Intubation Training Manikin (Adult) with Monitor RX-CPR184

### Description

- Pupil response: Bilateral pupils contrast observation;
- Carotid artery response: pinching pressure the ball to simulate carotid pulse
- Chest compressions
- Standard airway open
- Feasible three ways to open the airway: lift jaw with head back, lift the head and neck, lift jaw by hands
- Mouth resuscitation and ventilator-assisted breathing, effective artificial respiration has visible chest rise
- Simulate laryngeal mask airway and composite cannula can be inserted into the airway D-0009RX-CPR184 1,13,000.00 CPR Training
- Verify endotracheal intubation is correct or not: a catheter inserted into the trachea, lungs will swell (plastic bladder) when inflated; If the catheter is inserted by mistake oesophagus, stomach will swell when inflated (plastic pouch). Observed pulmonary respiration
- Artificial ventilation
- Oral or nasal intubation and oesophageal intubation
- Open the airway and practise from resuscitator bag to breathing mask or from resuscitator bag to the cannula
- Respiratory tract is simple and can be replaced
- Easy to carry, suitable for field training.



### PRODUCT RANGE

