

MODELSCART

Perineum Cutting and Suturing Training Simulator

RX-F3A



Description

- This model is designed for practice of perineum cutting and suturing.
- It simulates the anatomy of female vulva, including the Clitoris, Urethral Opening, Vaginal opening, Anus, Perineum, and Superficial muscles, and sphincter, etc.
- To make caring for and closing birth wounds as realistic as possible, the M. sphincter ani externus (outer anal sphincter muscle) can be palpated and sutured in the perineal area.
- Muscle and skin layers can be sutured separately.
- Can practice episiotomy: Median episiotomy, Mid-lateral episiotomy, Lateral episiotomy
- Can practice tear wounds: Vaginal tears, Labial tear
- The simulator is mounted on a plastic framework, easy disassembly.
- There are handles in the back and bottom of the simulator to adjust tension of the perineum

PRODUCT RANGE

