

MODELSCART

PICC LINE (Peripherally Insulated Central Catheter)

RX-L69C



Description

- The model is suitable for central venipuncture of internal carotid vein and sub-clavical vein and peripheral venipuncture (arm veins) of cephalic vein.
- The cervical veins and arm veins can be filled with simulative blood.
- Features:
 1. The model is a simulative right half of adult torso with right arm. There are apparent anatomical marks: clavicle, suprasternal notch, sternocleidomastoid muscle, pectoralis major muscle, ribs and deltoid muscle.
 2. Lifelike veins: superior vena cava, internal jugular vein, sub-clavical vein, cephalic vein, basilic vein and median cubital vein, etc.
 3. With apparent anatomical marks, the model can be used to train puncture and catheterization of internal carotid vein, sub-clavical vein, cephalic vein and heart catheterization.
 4. The skin and veins can be replaced. Students will feel an apparent "POP" when puncture needle has been thrust into the veins.
 5. The heart floating swan-ganz intubation can be exercised.

PRODUCT RANGE

