

MODELSCART

Infant CPR Training Model

XC-416



Description

- Clearly anatomic characteristic, realistic touch feeling and lifelike skin color, vivid appearance
- Can perform artificial respiratory, extracardial compression and airway open
- Operation mode: CPR training, can perform compression and inflation
- ▶ CPR Monitor functions:
 - Monitor shows inflation volume, correct inflation volume: 30ml- 50ml; Fast ventilation > 30ml within 0.27 seconds indicated by red light. When the inflation volume is excessive, proper, indicators will accordingly be yellow, green, and alarm will ring accordingly
 - Monitor shows compression depth: correct compression depth: approximately 13mm-20mm When the compression depth is proper or excessive, indicators will accordingly be
 - or excessive, indicators will accordingly be green, yellow, and alarm will ring accordingly
 - Monitor show compression position: Right compression position indicated by red light on the middle chest
 - Operation frequency: ≥ 100 times/ min
 - Ratio of compression and artificial respiratory: 30:2/single or 15:2/double
 - Operation cycle: after 2 valid inflation, then 5 cycles CPR operation according to the ratio of compression and inflation 30:2 or 15:2

PRODUCT RANGE

